






Tested disinfectants

The manufacturers of the following cleaning and disinfecting agents have carried out material suitability tests and confirmed compatibility with our brush strips. In some cases, the tests were carried out in high doses to simulate compatibility over a long period of time. All disinfectants should be dosed and used according to the manufacturer's recommendations.

Manufacturer		Tested disinfectants	Tested concentration in %	Duration of the test	Approval by manufacturer	short - term use	long - term use
	Biotec	BTS 2000	1,8	90 days			
		BTS 2001					
		Des AM blau					
		Virodex					
	Calvatis	Calgonit DS 683	5	19 days			
	Diversey	Suma San con. D10.1 conc	2	30 days			
		SaniPerfect VS61	2	70 days			
		Suredis VT1	3				
	Ecolab	P3 topax 66	5	28 days			
		P3 topax 990	3				
	Kersia	Deptal AS	2	33 days			
			20				
	Kesla	1+1 Wofasteril SC super	1	30 days			
	MENNO	MENNO® Clean	6	120 days			
		MENNO® Florades					
		A-QUASAN®	8				
	Stockmeier	Lerades Clp 601	2	30 days			
		Lerasept AC	2				
		Lerasept L420	2,5				
		Lerasept SOC	6				

Tested disinfectants

The manufacturers of the following cleaning and disinfecting agents have carried out material suitability tests and confirmed compatibility with our brush strips. In some cases, the tests were carried out in high doses to simulate compatibility over a long period of time. All disinfectants should be dosed and used according to the manufacturer's recommendations.

Manufacturer		Tested disinfectants	Tested concentration in %	Duration of the test	Approval by manufacturer	short-term use	long-term use
	Witty	Blaue Sandy	100	4 days			
		Goldene Sandy					
	ASiRAL	ASiRAL SRC	100	10 days	